



platts+Nisbett







www.plattsnisbett.com



Robert Nisbett, 1977.

Welcome to the Platts & Nisbett catalogue, a complete reference guide to our range of Surgical Instruments.

Our aim in producing this catalogue is that it should act as a comprehensive user guide, an invaluable resource when making your purchasing decision.

So, why choose Platts & Nisbett?

Customer Service - Our friendly staff are on hand to deal with your enquiries and support you every step of the way. We offer a prompt and reliable service, and will always do our best to meet your requirements.

Experience - We have been hand making Surgical Instruments since 1977, and have a wealth of manufacturing expertise at your disposal. We can produce bespoke instruments if our customers prefer something 'made to measure'.

Quality - According to leading surgeons the quality of our instruments is unsurpassed. Each one is manufactured to international quality standards and bears the CE Mark. Our ISO 13485 Management System ensures a consistent quality product is supplied first time, every time.

Reliability - All our instruments carry a Five Year Guarantee, meaning you can buy with confidence. Our instruments can last a lifetime if they are used and maintained with care, ensuring replacement costs are minimal.

Price - By dealing directly with the manufacturer you are able to benefit from significant cost savings, an important factor when trying to get the most out of your budget.

Thank you for taking the time to browse our catalogue. If you cannot find the product you require, please get in touch and we will be delighted to help you.


Alyson Nisbett
Director




Atraumatic Teeth




Debakey
1 x 2 Teeth



Debakey
2 x 3 Teeth



Cooley
2 x 2 Teeth



Dardik
2 x 3 Teeth



Longitudinal and Lateral Serrations



Longitudinal



Lateral


Measuring Guide

Please use the measuring guide below in order to assist you when making your selection.

To measure the overall length, please measure from the tip of the jaw to the bow as illustrated.

To measure the jaw length, please measure from the tip of the jaw to the bottom of the plain section as illustrated.

All measurements in this catalogue are in millimetres (mm).



Contents

Retractors

Pages 1-30



Clamps

Pages 31-138



Nasal Forceps

Pages 139-146



Bone Cutters

Pages 147-156



Orthopaedics

Pages 157-225



Index

Pages 227-236



Retractors



WEST RETRACTOR



PN0001



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0001	West Retractor	135	18	3 x 4 Sharp

WEISLANDER RETRACTOR



PN0014



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0014	Weislander Retractor	135	18	3 x 4 Blunt



PN0951



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN1094	Weitlaner Retractor	110	15	2 x 3 Sharp
PN0951	Weitlaner Retractor	110	15	2 x 3 Semi-Sharp
PN0952	Weitlaner Retractor	110	15	2 x 3 Blunt

TRIVERS RETRACTOR



PN0015



PN0016



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0015	Travers Retractor	210	26	4 x 5 Blunt
PN0016	Travers Retractor, Curved	200	26	4 x 5 Blunt



PN0012



PN0009



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0012	Mollison Retractor, Child	130	12	2 x 2 Sharp
PN0009	Mollison Retractor, Adult	170	15	4 x 4 Sharp



PN0007



1/1



PN0021



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0007	Norfolk & Norwich Retractor	220	50	4 x 5 Blunt
PN0021	Norfolk & Norwich Retractor, Heavy Rack	220	50	4 x 5 Blunt



PN0918



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0918	Norfolk & Norwich Retractor, Long Teeth	220	78	4 x 5 Blunt

JOLLS THYROID RETRACTOR



PN0212

CODE	DESCRIPTION	LENGTH
PN0212	Jolls Thyroid Retractor	155



PN0013



PN0989



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0013	Henderson Retractor	175	20	4 x 4 Sharp
PN0989	Henderson Retractor	175	20	4 x 4 Blunt



PN0008



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0008	Jefferson Retractor	130	16	3 x 4 Blunt

IRWIN RETRACTOR



PN0018



1/1


CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0018	Irwin Retractor	175	20	4 x 4 Sharp



PN0443



PN0444



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0443	Plester Jansen Retractor, Left Solid Blade	135	Left 28, Right 20	Solid x 3 Blunt
PN0444	Plester Jansen Retractor, Right Solid Blade	135	Left 20, Right 28	3 Blunt x Solid



PN0445



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0445	Plester Jansen Retractor	135	20	3 x 3 Blunt

WALLACE MODIFIED PLESTER JANSEN RETRACTOR



PN0410



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0410	Wallace Modified Plester Jansen Retractor	135	20	2 x 2 Semi-Sharp

CAWTHORNE RETRACTOR



PN0017



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0017	Cawthorne Retractor, Long Posterior Tooth	170	15	4 x 4 Sharp

MAWSON RETRACTOR



PN0019




1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0019	Mawson Retractor	140	12	3 x 3 Blunt



PN0557




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0557	Juniper Retractor	150	30	Blunt

JACKSON BURROWS RETRACTOR



PN0004




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0004	Jackson Burrows Retractor	180	20	Blunt



PN0005



1/1



PN0110




1/1


CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0005	Gelpi Retractor, Adult	185	20	Sharp
PN0110	Gelpi Spinal Retractor	280	78	Sharp




PN0067



PN0547



PN0067



PN0038

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0547	Gelpi Retractor	185	30	Sharp
PN0067	Gelpi Retractor	195	50	Sharp
PN0038	Gelpi Retractor	200	70	Sharp



PN0023



PN0125



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0023	Hardman Retractor, Equal Prongs	195	10	Sharp
PN0125	Hardman Retractor, Equal Prongs	215	50	Sharp



PN0022



PN0549



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0022	Hardman Retractor, Long Prong Left	215	Left 50, Right 30	Sharp
PN0549	Hardman Retractor, Long Prong Right	215	Left 30, Right 50	Sharp

CONE HEMI-LAMINECTOMY RETRACTOR



PN0447




PN0448



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0447	Cone Hemi-Laminectomy Retractor, Left	165	Left 17	Sharp
PN0448	Cone Hemi-Laminectomy Retractor, Right	165	Right 17	Sharp



PN0449



1/1



PN0450



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0449	Cone Hemi-Laminectomy Retractor, Left	255	Left 24	Sharp
PN0450	Cone Hemi-Laminectomy Retractor, Right	255	Right 24	Sharp

CONE RETRACTOR



PN0495



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0495	Cone Retractor	165	17	3 x 4 Sharp

CONE LAMINECTOMY RETRACTOR




PN0446



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0446	Cone Laminectomy Retractor	265	24	3 x 4 Sharp



PN0010



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0010	Harvey Jackson Laminectomy Retractor	305	45	3 x 3 Sharp

MAYO ADSON LAMINECTOMY RETRACTOR




PN0006




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0006	Mayo Adson Laminectomy Retractor	305	30	4 x 4 Blunt



PN0211



PN0440



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0211	Browse Retractor, Long Prong Left	195	Left 50, Right 25	4 x 4 Blunt
PN0440	Browse Retractor, Long Prong Right	195	Left 25, Right 50	4 x 4 Blunt



PN0441



1/1



PN0020



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0441	Talbot Orbital Retractor, Left	165	Left 40 & 52, Right 13	Interchangeable x 4 Sharp
PN0020	Talbot Orbital Retractor, Right	165	Left 13, Right 40 & 52	4 Sharp x Interchangeable

CAIRNS SCALP RETRACTOR



PN0115



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0115	Cairns Scalp Retractor	140	15	3 x 3 Sharp

R.N.T.N.E.H. RETRACTOR



PN0011



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0011	R.N.T.N.E.H. Retractor	170	15	3 x 4 Sharp



PN0324



1/1

Adult Blades



1/1

Sphincterotomy Blades



1/1

Centre Blade

CODE	DESCRIPTION	LENGTH
PN0324	Park Anal Retractor, with Adult, Sphincterotomy and Centre Blades	145



PN0558



PN0573

CODE	DESCRIPTION	LENGTH	BLADE DEPTH
PN0558	Eisenhammer Retractor	180	80
PN0573	Eisenhammer Retractor with Ratchet	180	80



PN0131

CODE	DESCRIPTION	LENGTH	BLADE DEPTH
PN0131	Farrands Retractor	180	140



PN0041

CODE	DESCRIPTION
PN0041	Sellors Retractor



PN1058



PN1059

CODE	DESCRIPTION	LENGTH
PN1058	Harefield Sternal Blades	78 x 30
PN1059	Harefield Sternal Blades	120 x 30

IMA RETRACTOR



PN0305

CODE	DESCRIPTION
PN0305	IMA Retractor

PHLEBECTOMY HOOK



PN0801

CODE	DESCRIPTION
PN0801	Phlebectomy Hook No. 1
PN0802	Phlebectomy Hook No. 2
PN0803	Phlebectomy Hook No. 3

Clamps



COOLEY CROSS ACTION BULLDOG CLAMP




PN1085




PN0118



PN0592



PN1084





PN1085

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0118	Cooley Cross Action Bulldog Clamp, Straight	80	24	Cooley
PN0592	Cooley Cross Action Bulldog Clamp, Straight	90	35	Cooley
PN1084	Cooley Cross Action Bulldog Clamp, Straight	105	50	Cooley
PN1085	Cooley Cross Action Bulldog Clamp, Straight	124	69	Cooley



PN1087



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0310	Cooley Cross Action Bulldog Clamp, Curved	77	20	Cooley
PN0594	Cooley Cross Action Bulldog Clamp, Curved	86	30	Cooley
PN1086	Cooley Cross Action Bulldog Clamp, Curved	98	44	Cooley
PN1087	Cooley Cross Action Bulldog Clamp, Curved	115	63	Cooley

DEBAKEY CROSS ACTION BULLDOG CLAMP



PN0591



1/1

PN0545



1/1

PN0309



1/1

PN1093



1/1


PN0590



1/1

PN0591

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0545	Debakey Cross Action Bulldog Clamp, Straight	80	24	Debakey 1 x 2
PN0309	Debakey Cross Action Bulldog Clamp, Straight	90	34	Debakey 1 x 2
PN1093	Debakey Cross Action Bulldog Clamp, Straight	100	44	Debakey 1 x 2
PN0590	Debakey Cross Action Bulldog Clamp, Straight	105	49	Debakey 1 x 2
PN0591	Debakey Cross Action Bulldog Clamp, Straight	125	69	Debakey 1 x 2

**PN0599****PN0308****PN0546****PN0554****PN0599****PN0057**

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0308	Debakey Cross Action Bulldog Clamp, Curved	75	22	Debakey 1 x 2
PN0546	Debakey Cross Action Bulldog Clamp, Curved	85	30	Debakey 1 x 2
PN0554	Debakey Cross Action Bulldog Clamp, Curved	95	41	Debakey 1 x 2
PN0599	Debakey Cross Action Bulldog Clamp, Curved	100	44	Debakey 1 x 2
PN0057	Debakey Cross Action Bulldog Clamp, Curved	115	63	Debakey 1 x 2

COOLEY PAEDIATRIC OCCLUSION CLAMP




PN0182



1/1

PN0181



1/1

PN0183




1/1

PN0182



1/1

PN0184



1/1

PN0185

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0181	Cooley Paediatric Occlusion Clamp, Straight	143	37	Cooley
PN0183	Cooley Paediatric Occlusion Clamp, Angled Jaw	143	37	Cooley
PN0182	Cooley Paediatric Occlusion Clamp, Angled Shank	144	37	Cooley
PN0184	Cooley Paediatric Occlusion Clamp, Spoon Jaw	143	35	Cooley
PN0185	Cooley Paediatric Occlusion Clamp, Satinsky Jaw	140	35	Cooley



PN0550



PN1092 1/1



PN0550 1/1



PN0551 1/1




PN0134 1/1


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1092	Debakey Paediatric Occlusion Clamp, Straight	156	38	Debakey 1 x 2
PN0550	Debakey Paediatric Occlusion Clamp, Angled Jaw	155	36	Debakey 1 x 2
PN0551	Debakey Paediatric Occlusion Clamp, Angled Shank	155	37	Debakey 1 x 2
PN0134	Debakey Paediatric Occlusion Clamp, Spoon Jaw	155	34	Debakey 1 x 2




PN0318



PN0272



PN0336



PN0318

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0272	Cooley Derra Clamp, Small	170	46	Cooley
PN0336	Cooley Derra Clamp, Medium	175	50	Cooley
PN0318	Cooley Derra Clamp, Large	175	53	Cooley



PN0031




1/1

PN0029



1/1

PN0030



1/1

PN0031

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0029	Debakey Derra Clamp, Small	165	42	Debakey 1 x 2
PN0030	Debakey Derra Clamp, Medium	171	48	Debakey 1 x 2
PN0031	Debakey Derra Clamp, Large	173	50	Debakey 1 x 2

COOLEY BULLDOG CLAMP, RING HANDLE




PN0342




1/1

PN0176



1/1

PN0341



1/1

PN0178



1/1

PN0179



1/1

PN0342

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0176	Cooley Bulldog Clamp, Ring Handle, Straight	130	48	Cooley
PN0341	Cooley Bulldog Clamp, Ring Handle, Angled	127	46	Cooley
PN0178	Cooley Bulldog Clamp, Ring Handle, Curved	115	36	Cooley
PN0179	Cooley Bulldog Clamp, Ring Handle, Doubled Angled	127	46	Cooley
PN0342	Cooley Bulldog Clamp, Ring Handle, Satinsky Jaws	120	35	Cooley




PN0375



1/1

PN0376




1/1

PN0188



1/1

PN0189



1/1

PN0569



1/1

PN0375

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0376	Debakey Bulldog Clamp, Ring Handle, Straight	135	50	Debakey 1 x 2
PN0188	Debakey Bulldog Clamp, Ring Handle, 45°	130	46	Debakey 1 x 2
PN0189	Debakey Bulldog Clamp, Ring Handle, 90°	105	38	Debakey 1 x 2
PN0569	Debakey Bulldog Clamp, Ring Handle, Curved Jaws	130	46	Debakey 1 x 2
PN0375	Debakey Bulldog Clamp, Ring Handle, Double Curved	130	47	Debakey 1 x 2

WOLFE MODIFIED COOLEY BULLDOG CLAMP



PN0143

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0143	Wolfe Modified Cooley Bulldog Clamp, Ring Handle, Double Angled	160	44	Cooley

COOLEY VESSEL CLAMP




PN0144


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0144	Cooley Vessel Clamp	165	50	Debakey 1 x 2



PN0370



PN0026



PN0027



PN0370

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0026	Debakey Miniature Clamp, 30°	165	47	Debakey 1 x 2
PN0027	Debakey Miniature Clamp, 60°	160	44	Debakey 1 x 2
PN0370	Debakey Miniature Clamp, 90°	146	41	Debakey 1 x 2


COOLEY MULTI PURPOSE CLAMP




PN0035



PN0034



PN0035



PN0036

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0034	Cooley Multi Purpose Clamp, 30°	165	47	Cooley
PN0035	Cooley Multi Purpose Clamp, 60°	160	43	Cooley
PN0036	Cooley Multi Purpose Clamp, 90°	148	38	Cooley




PN0337



1/1



PN0371



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0337	Debakey Multi Purpose Clamp No. 1, 60°	245	68	Debakey 1 x 2
PN0371	Debakey Multi Purpose Clamp No. 2, 60°	205	60	Debakey 1 x 2

DEBAKEY MULTI PURPOSE CLAMP




PN0154



1/1



PN0045



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0154	Debakey Multi Purpose Clamp No. 3, 60°	242	68	Debakey 1 x 2
PN0045	Debakey Multi Purpose Clamp No. 4, 60°	240	75	Debakey 1 x 2




PN0372



1/1



PN0373



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0372	Debakey Multi Purpose Clamp No. 5, 60°	280	87	Debakey 1 x 2
PN0373	Debakey Multi Purpose Clamp No. 6, 60°	310	125	Debakey 1 x 2

DEBAKEY MULTI PURPOSE CLAMP



PN0540



PN0083



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0540	Debakey Multi Purpose Clamp No. 1, 90°	230	60	Debakey 1 x 2
PN0083	Debakey Multi Purpose Clamp No. 2, 90°	187	57	Debakey 1 x 2



PN0355




PN0356



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0355	Debakey Multi Purpose Clamp No. 3, 90°	220	70	Debakey 1 x 2
PN0356	Debakey Multi Purpose Clamp No. 4, 90°	270	110	Debakey 1 x 2

DARDIK MULTI PURPOSE CLAMP




PN0122




1/1

PN0351



1/1

PN0349



1/1

PN0122



1/1

PN0972



1/1

PN0124



1/1

PN0350

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0351	Dardik Multi Purpose Clamp, Small, Straight	150	35	Dardik
PN0349	Dardik Multi Purpose Clamp, Small, 15°	149	34	Dardik
PN0122	Dardik Multi Purpose Clamp, Small, 30°	148	33	Dardik
PN0972	Dardik Multi Purpose Clamp, Small, 45°	145	32	Dardik
PN0124	Dardik Multi Purpose Clamp, Small, 50°	142	31	Dardik
PN0350	Dardik Multi Purpose Clamp, Small, 90°	130	28	Dardik




PN0353




1/1

PN0492




1/1

PN0352



1/1

PN0116



1/1

PN0117



1/1

PN0353



1/1


PN0119

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0492	Dardik Multi Purpose Clamp, Large, Straight	170	55	Dardik
PN0352	Dardik Multi Purpose Clamp, Large, 15°	169	54	Dardik
PN0116	Dardik Multi Purpose Clamp, Large, 30°	165	53	Dardik
PN0117	Dardik Multi Purpose Clamp, Large, 45°	161	51	Dardik
PN0353	Dardik Multi Purpose Clamp, Large, 50°	158	50	Dardik
PN0119	Dardik Multi Purpose Clamp, Large, 90°	140	45	Dardik

CASTENADA CLAMP



PN1072




CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1069	Castenada Clamp	110	35	Debakey 1 x 2
PN1070	Castenada Clamp	117	37	Debakey 1 x 2
PN0598	Castenada Clamp	126	42	Debakey 1 x 2
PN1071	Castenada Clamp	110	25	Debakey 1 x 2
PN1072	Castenada Clamp	120	37	Debakey 1 x 2
PN0177	Castenada Clamp	120	35	Debakey 1 x 2
PN1073	Castenada Clamp	120	35	Debakey 1 x 2
PN1074	Castenada Clamp	125	43	Debakey 1 x 2
PN1075	Castenada Clamp	120	40	Debakey 1 x 2
PN1076	Castenada Clamp	115	30	Debakey 1 x 2




PN0063




PN0059




PN0494



PN1077



PN0062




PN0063

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0059	Neonatal Miniature Clamp, Straight	128	47	Lateral
PN0494	Neonatal Miniature Clamp, 30°	127	46	Lateral
PN1077	Neonatal Miniature Clamp, 50°	125	43	Lateral
PN0062	Neonatal Miniature Clamp, 90°	112	36	Lateral
PN0063	Neonatal Miniature Clamp, Spoon	123	42	Lateral

DEBAKEY PERIPHERAL CLAMP




PN0025




1/1



PN0033



1/1

**PN0382****1/1**

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0025	Debakey Peripheral Clamp, Angled	183	70	Debakey 1 x 2
PN0033	Debakey Peripheral Clamp, S Shape	200	83	Debakey 1 x 2
PN0382	Debakey Peripheral Vascular Clamp, Deep Curve	160	47	Debakey 1 x 2

LELAND JONES PERIPHERAL VASCULAR CLAMP




PN0555




1/1

PN0064




1/1

PN0128



1/1

PN1096



1/1

PN0555

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0064	Leland Jones Peripheral Vascular Clamp, Straight	195	79	Debakey 1 x 2
PN0128	Leland Jones Peripheral Vascular Clamp, 15°	195	77	Debakey 1 x 2
PN1096	Leland Jones Peripheral Vascular Clamp, 30°	190	75	Debakey 1 x 2
PN0555	Leland Jones Peripheral Vascular Clamp, 45°	186	75	Debakey 1 x 2



PN0587




CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0587	Debakey Renal Clamp	198	80	Debakey 1 x 2

DEBAKEY ILIAC CLAMP



PN0294





CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0294	Debakey Iliac Clamp	200	67	Debakey 1 x 2




PN1080

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1078	Debakey Tangential Clamp No. 1	195	46	Debakey 1 x 2
PN1079	Debakey Tangential Clamp No. 2	255	55	Debakey 1 x 2
PN1080	Debakey Tangential Clamp No. 3	260	65	Debakey 1 x 2
PN1081	Debakey Tangential Clamp No. 4	270	75	Debakey 1 x 2
PN1082	Debakey Tangential Clamp No. 5	280	90	Debakey 1 x 2





PN0101



1/1



PN0384



1/1



PN0099



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0101	Debakey Aneurysm Clamp, Shallow Curve	260	70	Debakey 1 x 2
PN0384	Debakey Aneurysm Clamp, Medium Curve	255	68	Debakey 1 x 2
PN0099	Debakey Aneurysm Clamp, Deep Curve	245	59	Debakey 1 x 2

DEBAKEY AORTIC ANEURYSM CLAMP




PN0358



1/1




PN0102



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0358	Debakey Aortic Aneurysm Clamp	235	84	Debakey 1 x 2
PN0102	Debakey Aortic Aneurysm Clamp	270	93	Debakey 1 x 2



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0107	Debakey Aortic Aneurysm Clamp	310	113	Debakey 1 x 2

DEBAKEY AORTIC ANEURYSM CLAMP




PN0103

1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0103	Debakey Aortic Aneurysm Clamp	310	113	Debakey 2 x 3



PN0024




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0024	Debakey Aortic Occlusion Clamp	195	80	Debakey 1 x 2

MATTOX AORTA CLAMP



PN0180



1/1


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0180	Mattox Aorta Clamp	258	90	Debakey 1 x 2



PN0085



1/1



PN0106



1/1



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0085	Debakey Aortic Clamp, Full Curve	270	80	Debakey 1 x 2
PN0106	Debakey Aortic Clamp, Spoon	270	80	Debakey 2 x 3
PN0105	Debakey Aortic Clamp, S Shape	270	78	Debakey 2 x 3

MORRIS AORTA CLAMP



PN0301



1/1



PN0407



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0301	Morris Aorta Clamp	210	85	Debakey 2 x 3
PN0407	Morris Aorta Clamp	282	87	Debakey 2 x 3



PN0164



1/1



PN0746



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0164	Lambert Kay Clamp	233	49	Debakey 1 x 2
PN0746	Lambert Kay Clamp	253	50	Debakey 1 x 2


BECK AORTIC CLAMP



PN0068



PN0068



PN0710

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0068	Beck Aortic Clamp	195	42	Debaeye 1 x 2
PN0710	Beck Aortic Clamp	205	60	Debaeye 1 x 2



PN0140



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0140	Satinsky Clamp, Baby	165	47	Debakey 1 x 2

SATINSKY CLAMP



PN0071



PN0069 1/1

PN0070 1/1

PN0071 1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0069	Satinsky Clamp, Small	235	50	Debakey 1 x 2
PN0070	Satinsky Clamp, Medium	240	55	Debakey 1 x 2
PN0071	Satinsky Clamp, Large	250	65	Debakey 1 x 2



PN0166

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0166	Ross Vena Cava Clamp, 26 Charriere	215	72	Cooley
PN0647	Ross Vena Cava Clamp, 28 Charriere	215	72	Cooley
PN0584	Ross Vena Cava Clamp, 32 Charriere	215	72	Cooley

JAVID CAROTID ARTERY CLAMP



PN0532



PN0291

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0532	Javid Carotid Artery Clamp, Small	170	55
PN0291	Javid Carotid Artery Clamp, Large	185	75




PN1088



1/1



PN1090



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1088	Debaquey Bainbridge Clamp, Straight	165	46	Debaquey 1 x 2
PN1090	Debaquey Bainbridge Clamp, Curved	165	45	Debaquey 1 x 2

DEBAKEY BAINBRIDGE CLAMP



PN1089



1/1



PN0359



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1089	Debakey Bainbridge Clamp, Straight	185	65	Debakey 1 x 2
PN0359	Debakey Bainbridge Clamp, Curved	185	68	Debakey 1 x 2



PN0361



1/1



PN0360



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0361	Debakey Coarctation Clamp, Straight	170	55	Debakey 1 x 2
PN0360	Debakey Coarctation Clamp, Angled Shank	170	55	Debakey 1 x 2

DEBAKEY COARCTATION CLAMP



PN0077



1/1



PN0078



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0077	Debakey Coarctation Clamp, Straight	220	67	Debakey 1 x 2
PN0078	Debakey Coarctation Clamp, Angled Shank	215	67	Debakey 1 x 2



PN0507

1/1



PN0163


1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0507	Debakey Coarctation Clamp, Straight	257	80	Debakey 1 x 2
PN0163	Debakey Coarctation Clamp, Angled Shank	255	80	Debakey 1 x 2

DEBAKEY COARCTATION CLAMP



PN0052



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0052	Debakey Coarctation Clamp, Curved	165	50	Debakey 1 x 2




PN0705



PN0571



PN0039



PN0097



PN0705

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0571	Potts Coarctation Clamp, Straight	209	59	Lateral
PN0039	Potts Coarctation Clamp, Angled Shank	205	58	Lateral
PN0097	Potts Coarctation Clamp, Curved	205	58	Lateral
PN0705	Potts Coarctation Clamp, Spoon	203	62	Lateral


CRAFOORD COARCTATION CLAMP




PN0347



1/1



PN0157




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0347	Crafoord Coarctation Clamp, Straight	255	78	Longitudinal
PN0157	Crafoord Coarctation Clamp, Slight Curve	253	75	Longitudinal




PN0348



1/1



PN0295




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0348	Crafoord Dissecting Clamp, Medium Curve	245	68	Lateral
PN0295	Crafoord Dissecting Clamp, Full Curve	260	85	Lateral

GLOVER CLAMP



PN0079



1/1



PN0393




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0079	Glover Clamp, Curved	210	60	Debakey 1 x 2
PN0393	Glover Clamp, Spoon	215	60	Debakey 1 x 2



PN0096



1/1



PN0394




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0096	Glover Clamp, Curved	250	73	Debakey 1 x 2
PN0394	Glover Clamp, Spoon	257	77	Debakey 1 x 2

BROCK AURICULAR CLAMP




PN0088



1/1
PN0340



1/1
PN0088




1/1
PN0339

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0340	Brock Auricular Clamp, Shallow Curve	232	82	Brock
PN0088	Brock Auricular Clamp, Medium Curve	230	80	Brock
PN0339	Brock Auricular Clamp, Deep Curve	213	70	Brock



PN0402



1/1



PN0401



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0402	Lees Bronchus Clamp, Slight Angled Jaw	245	65	Debakey 1 x 2
PN0401	Lees Bronchus Clamp, 90°	223	57	Debakey 1 x 2




PN0273



1/1
PN0274



1/1
PN0273




1/1
PN0243

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0274	Price Thomas Bronchus Clamp, Plain	220	50	Longitudinal
PN0273	Price Thomas Bronchus Clamp, Teeth on Left	220	50	Longitudinal
PN0243	Price Thomas Bronchus Clamp, Teeth on Right	220	50	Longitudinal




PN0086



1/1



PN0293




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0086	Semb Ligature Passing Forceps, Child, Full Curve	210	55
PN0293	Semb Ligature Passing Forceps, Adult, Full Curve	245	77

SEMB LIGATURE PASSING FORCEPS



PN0418



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0418	Semb Ligature Passing Forceps, Adult, Slight Curve	245	73



PN0417



1/1



PN0046



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0417	Debakey Modified Semb Ligature Passing Forceps, Full Curve, Serrated Tips	264	68	Debakey 1 x 2
PN0046	Debakey Ligature Forceps, Full Curve	170	50	Debakey 1 x 2



PN0426



1/1



PN0425



1/1



PN0424


1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0426	Ulrich Clamp, Small	260	64	Lateral
PN0425	Ulrich Clamp, Medium	270	75	Lateral
PN0424	Ulrich Clamp, Large	280	87	Lateral

SWINNEY CLAMP




PN0420



1/1



PN0419



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0420	Swinney Clamp, Child	150	45	Debakey 1 x 2
PN0419	Swinney Clamp, Adult	235	60	Debakey 1 x 2



PN0267

1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0267	Tubing Clamp, Heavy Pattern	200	48	Lateral

TUBING CLAMP



PN0160



1/1

PN0159



1/1

PN0160



1/1

PN0161

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0159	Tubing Clamp, Wedge Jaw	150	34
PN0160	Tubing Clamp, Wedge Jaw	182	38
PN0161	Tubing Clamp, Wedge Jaw	207	42



PN0422



1/1

PN0421



1/1

PN0422




1/1

PN0423

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0421	Tubing Clamp, Wedge Jaw, With Safety Guard	150	34
PN0422	Tubing Clamp, Wedge Jaw, With Safety Guard	182	38
PN0423	Tubing Clamp, Wedge Jaw, With Safety Guard	207	42



PN0436



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0436	Lahey Artery Forceps	200	40	Lateral



PN0175



1/1



PN0174




1/1


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0175	Lahey Artery Forceps	200	40	Longitudinal
PN0174	Lahey Artery Forceps	235	54	Longitudinal



PN0322



PN0403



PN0322

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0403	Lloyd Davis Artery Forceps, Slight Curve	257	53	Lateral
PN0322	Lloyd Davis Artery Forceps, Deep Curve	255	50	Lateral



1/1

PN0973

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0973	Nelson Roberts Forceps	235	55	Lateral

SAWTELL FORCEPS



1/1

PN0965

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0965	Sawtell Forcpes	230	40	Lateral

ROBERTS ARTERY FORCEPS



PN0191

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0191	Roberts Artery Forceps	300	55	Lateral

MOYNIHAN CLAMP




PN0152


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0152	Moynihan Clamp, Curved	225	40	Lateral



PN0098



PN0100




CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0098	Waterston Dissecting Forceps, Angled	155	35	Lateral
PN0100	Waterston Dissecting Forceps, Curved	167	41	Lateral


WATERSTON DISSECTING FORCEPS



PN0711



PN0696



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0711	Waterston Dissecting Forceps, Angled, Delicate Model	155	33	Lateral
PN0696	Waterston Dissecting Forceps, Curved, Delicate Model	155	40	Lateral



1/1

PN0706



1/1


PN0694

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0706	Waterston Dissecting Forceps, Angled, Delicate Model	170	43	Lateral
PN0694	Waterston Dissecting Forceps, Curved, Delicate Model	170	41	Lateral


WATERSTON DISSECTING FORCEPS



PN0949




PN0948




CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0949	Waterston Dissecting Forceps, Angled, with Bulbous Tip	155	35	Lateral
PN0948	Waterston Dissecting Forceps, Curved, with Bulbous Tip	165	45	Lateral



PN0585



PN0578




CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0585	Waterston Dissecting Forceps, Angled	203	58	Lateral
PN0578	Waterston Dissecting Forceps, Curved	203	58	Lateral


WATERSTON VESSEL CLAMP



PN0583



1/1



PN0302




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0583	Waterston Vessel Clamp, Delicate Model	165	47	Cooley
PN0302	Waterston Vessel Clamp, Delicate Model	180	57	Cooley



PN0104



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0104	Lincoln Dissecting Forceps	115	15	Lateral

COLLINGWOOD STEWART HERNIA FORCEPS



PN0095

CODE	DESCRIPTION	LENGTH
PN0095	Collingwood Stewart Hernia Forceps	125

KING VEIN GRAFT PRESS FORCEPS



PN0541

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0541	King Vein Graft Press Forceps	145	20	Circular

PURSE STRING FORCEPS




PN0087

CODE	DESCRIPTION	LENGTH
PN0087	Purse String Forceps	240



PN0335




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0335	Herrick Kidney Pedicle Clamp	260	82	Debakey 1 x 2

FISTULA CLAMP



PN0383



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0383	Fistula Clamp	117	35	Debakey 1 x 2

DOYEN CLAMP



PN0628



1/1



PN0306




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0628	Doyen Clamp, Baby, Straight	177	64	Longitudinal
PN0306	Doyen Clamp, Baby, Curved	175	60	Longitudinal




PN0627



1/1




PN0679



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0627	Doyen Clamp, Adult, Straight	230	100	Longitudinal
PN0679	Doyen Clamp, Adult, Curved	230	98	Longitudinal


DOYEN CLAMP




PN0387



1/1




PN0129



1/1


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0387	Doyen Clamp, Baby, Straight	177	62	Debakey 2 x 3
PN0129	Doyen Clamp, Baby, Curved	175	60	Debakey 2 x 3




PN0388



1/1



PN0130



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0388	Doyen Clamp, Adult, Straight	230	102	Debakey 2 x 3
PN0130	Doyen Clamp, Adult, Curved	228	98	Debakey 2 x 3



PN0552



1/1

PN0552




1/1

PN0300


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0552	Parker Kerr Clamp, Small, Straight	242	65	Longitudinal
PN0300	Parker Kerr Clamp, Small, Curved	238	61	Longitudinal



PN0270



PN0299



PN0270

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0299	Parker Kerr Clamp, Large, Straight	258	83	Longitudinal
PN0270	Parker Kerr Clamp, Large, Curved	255	80	Longitudinal



PN0213



PN0414

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0213	Payr Intestinal Clamp	202	58	Longitudinal
PN0414	Payr Intestinal Clamp	250	106	Longitudinal



PN0123



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0123	Tanner (Lane) Twin Clamp, Straight	313	165	Longitudinal



PN0400



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0400	Lane Twin Clamp, Straight	307	157	Longitudinal




PN0399

1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0399	Lane Twin Clamp, Curved	305	155	Longitudinal



PN1095




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN1095	Geary Grant Forceps, Small	220	43	Longitudinal



PN0307



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0307	Geary Grant Forceps, Large	225	50	Longitudinal

DEBAKEY HYSTERECTOMY CLAMP




PN0142



1/1




PN0365



1/1


CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0142	Debakey Hysterectomy Clamp, Straight	220	54	Debakey 2 x 3
PN0365	Debakey Hysterectomy Clamp, Slight Curved Jaws	220	54	Debakey 2 x 3




PN0366



1/1




PN0364



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0366	Debakey Hysterectomy Clamp, Deep Curved Jaws	215	49	Debakey 2 x 3
PN0364	Debakey Hysterectomy Clamp, Angled Jaws	205	44	Debakey 2 x 3


DEBAKEY HYSTERECTOMY CLAMP




PN0146



1/1



PN0368




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0146	Debakey Hysterectomy Clamp, Straight	243	55	Debakey 2 x 3
PN0368	Debakey Hysterectomy Clamp, Slight Curved Jaws	243	54	Debakey 2 x 3



PN0369



1/1



PN0367




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0369	Debakey Hysterectomy Clamp, Deep Curved Jaws	240	50	Debakey 2 x 3
PN0367	Debakey Hysterectomy Clamp, Angled Jaws	225	43	Debakey 2 x 3


MAINGOT HYSTERECTOMY FORCEPS



PN0405



PN0405



PN0153

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0405	Maingot Hysterectomy Forceps, Small, Straight, Tenaculum Point	210	59
PN0153	Maingot Hysterectomy Forceps, Small, Curved, Tenaculum Point	205	57



PN0406



1/1

PN0319



1/1

PN0406

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0319	Maingot Hysterectomy Forceps, Large, Straight, Tenaculum Point	251	61
PN0406	Maingot Hysterectomy Forceps, Large, Curved, Tenaculum Point	248	59

GWILLIAMS HYSTERECTOMY CLAMP



PN0139



1/1


PN0395



1/1

PN0139

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0395	Gwilliams Hysterectomy Clamp, Straight, 1 x 2 Tenaculum Point	208	57	Longitudinal
PN0139	Gwilliams Hysterectomy Clamp, Curved, 1 x 2 Tenaculum Point	205	55	Longitudinal




PN0032



1/1



PN0138




1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0032	Rogers Hysterectomy Clamp, Straight	215	50	Debakey 2 x 3
PN0138	Rogers Hysterectomy Clamp, Curved	214	49	Debakey 2 x 3

READ STALLWORTHY HYSTERECTOMY CLAMP




PN0141



1/1



PN0084



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0141	Read Stallworthy Hysterectomy Clamp (Chelsea Pattern), Tenaculum Point	205	42	Lateral
PN0084	Read Stallworthy Hysterectomy Clamp, Tenaculum Point	224	47	Lateral



PN0564



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0564	Howkin Hysterectomy Clamp, 1 x 2 Tenaculum Point	212	46	Longitudinal

TEALE VULSELLUM FORCEPS



PN0357

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0357	Teale Vulsellum Forceps	225	55	3 x 4

FINCH EXCISION CLAMP



PN0390

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0390	Finch Excision Clamp	235	95	Hollow

LLOYD DAVIS RECTAL OCCLUSION FORCEPS



PN0061

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0061	Lloyd Davis Rectal Occlusion Forceps	330	95	Longitudinal



PN0396



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0396	Houghton Excision Clamp	255	85	Longitudinal

EINSELL FOSS COLON CLAMP



PN0389



1/1

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0389	Einsell Foss Colon Clamp	310	90	Debakey 1 x 2

HAYES COLON RESECTION CLAMP



PN0108



CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0108	Hayes Colon Resection Clamp	270	80	Debakey 1 x 2
PN0109	Hayes Colon Resection Clamp	285	100	Debakey 1 x 2



PN0028



1/1

PN0028




1/1


PN0055

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0028	Goldberg Anterior Resection Clamp	335	75	Debakey 1 x 2
PN0055	Goldberg Anterior Resection Clamp	340	85	Debakey 1 x 2

McDOUGALL PROSTATECTOMY CLAMP



PN0298



PN0297

CODE	DESCRIPTION	LENGTH	JAW LENGTH	TEETH
PN0298	McDougall Prostatectomy Clamp, Left	255	40	Lateral
PN0297	McDougall Prostatectomy Clamp, Right	255	40	Lateral


Nasal Forceps





PN0428

CODE	DESCRIPTION	LENGTH	BLADE LENGTH
PN0427	Killian Nasal Speculum	135	50
PN0428	Killian Nasal Speculum	135	55
PN0429	Killian Nasal Speculum	135	65
PN0430	Killian Nasal Speculum	135	75
PN0431	Killian Nasal Speculum	135	90





PN0040

CODE	DESCRIPTION	LENGTH	BLADE LENGTH
PN0040	Cottle Nasal Speculum - Satin	140	20
PN0092	Cottle Nasal Speculum - Satin	140	32
PN0290	Cottle Nasal Speculum - Satin	140	41
PN0320	Cottle Nasal Speculum - Satin	140	60
PN1136	Cottle Nasal Speculum - Black	140	20
PN1137	Cottle Nasal Speculum - Black	140	32
PN1138	Cottle Nasal Speculum - Black	140	41
PN1139	Cottle Nasal Speculum - Black	140	60



1/1
PN0040



1/1
PN0092



1/1
PN0290



1/1
PN0320



PN0167



1/1
PN0167



1/1
PN0168



1/1
PN0169

CODE	DESCRIPTION	SHAPE	LENGTH	JAW WIDTH
PN0167	Luc Nasal Forceps	Oval	195	6
PN0168	Luc Nasal Forceps	Oval	195	7.5
PN0169	Luc Nasal Forceps	Oval	195	9

**PN0170****PN0170****PN0171****PN0172**

CODE	DESCRIPTION	SHAPE	LENGTH	JAW WIDTH
PN0170	Luc Nasal Forceps	Heart	195	6
PN0171	Luc Nasal Forceps	Heart	195	7.5
PN0172	Luc Nasal Forceps	Heart	195	9



PN0268



1/1

PN0268

CODE	DESCRIPTION	LENGTH	JAW WIDTH
PN0268	Irwin Moore Nasal Forceps	200	4

Bone Cutters



STAMM BONE CUTTING FORCEPS



PN0279

CODE	DESCRIPTION	LENGTH
PN0279	Stamm Bone Cutting Forceps	160

PATON BONE CUTTING FORCEPS



PN0432

CODE	DESCRIPTION	LENGTH
PN0432	Paton Bone Cutting Forceps	140

**PN0278****PN0338**

CODE	DESCRIPTION	LENGTH
PN0278	McIndoe Bone Cutting Forceps, Compound Action	170
PN0338	McIndoe Bone Cutting Forceps, Compound Action	180

LISTON BONE CUTTING FORCEPS



PN0332



PN0236

CODE	DESCRIPTION	LENGTH
PN0332	Liston Bone Cutting Forceps, Straight Jaw	175
PN0236	Liston Bone Cutting Forceps, Straight Jaw	200



PN0094



PN0437

CODE	DESCRIPTION	LENGTH
PN0094	Liston Bone Cutting Forceps, Angled Jaw	165
PN0437	Liston Bone Cutting Forceps, Angled Jaw	200

GLASGOW RONGEUR



PN0275

CODE	DESCRIPTION	LENGTH	BITE
PN0275	Glasgow Rongeur	165	4

PENNYBACKER RONGEUR



PN0333

CODE	DESCRIPTION	LENGTH	BITE
PN0333	Pennybacker Rongeur	205	4



PN0589

CODE	DESCRIPTION	LENGTH
PN0589	Ward Bone Cutting Forceps	170

WILMS RONGEUR



PN0619

CODE	DESCRIPTION	LENGTH	BITE
PN0619	Wilms Rongeur	205	4

LUER JANSEN RONGEUR



PN0276

CODE	DESCRIPTION	LENGTH	BITE
PN0276	Luer Jansen Rongeur, Compound Action	195	4

CAWTHORNE BONE NIBBLING FORCEPS



PN0544

CODE	DESCRIPTION	LENGTH	BITE
PN0544	Cawthorne Bone Nibbling Forceps	160	3



PN0237

CODE	DESCRIPTION	LENGTH	BITE
PN0237	Jansen Middleton Septum Forceps, Compound Action	185	4

OLIVECRONA RONGEUR



PN0879

CODE	DESCRIPTION	LENGTH	BITE
PN0879	Olivecrona Rongeur, Angled to Side, Compound Action	230	6

HORSLEY LISTON BONE CUTTING FORCEPS



PN0238

CODE	DESCRIPTION	LENGTH
PN0238	Horsley Liston Bone Cutting Forceps, Compound Action	280

VAN HAVRE RONGEUR



PN0435

CODE	DESCRIPTION	LENGTH	BITE
PN0435	Van Havre Rongeur, Compound Action	270	10

Orthopaedics




 A straight chisel with a wide, flat blade and a smooth, light-colored wooden handle. The blade has a small logo near the handle.

PN0463


 A curved chisel with a narrow, slightly curved blade and a smooth, light-colored wooden handle. The blade has a small logo near the handle.

PN0137

CODE	DESCRIPTION	LENGTH	WIDTH
PN0460	Capener Chisel, Straight	260	6
PN0977	Capener Chisel, Straight	260	8
PN0461	Capener Chisel, Straight	260	10
PN0462	Capener Chisel, Straight	260	12
PN0463	Capener Chisel, Straight	260	15
PN0464	Capener Chisel, Straight	260	20
PN0465	Capener Chisel, Straight	260	25
PN0978	Capener Chisel, Straight	260	30
PN0137	Capener Chisel, Curved	260	6
PN0979	Capener Chisel, Curved	260	8
PN0980	Capener Chisel, Curved	260	10
PN0981	Capener Chisel, Curved	260	12
PN0982	Capener Chisel, Curved	260	15
PN0983	Capener Chisel, Curved	260	20
PN0984	Capener Chisel, Curved	260	25
PN0985	Capener Chisel, Curved	260	30

PN0451



PN0458



CODE	DESCRIPTION	LENGTH	WIDTH
PN0451	Capener Osteotome, Straight	260	6
PN0089	Capener Osteotome, Straight	260	8
PN0313	Capener Osteotome, Straight	260	10
PN0314	Capener Osteotome, Straight	260	12
PN0452	Capener Osteotome, Straight	260	15
PN0315	Capener Osteotome, Straight	260	20
PN0311	Capener Osteotome, Straight	260	25
PN0312	Capener Osteotome, Straight	260	30
PN0702	Capener Osteotome, Straight	260	40
PN0453	Capener Osteotome, Curved	260	6
PN0987	Capener Osteotome, Curved	260	8
PN0454	Capener Osteotome, Curved	260	10
PN0455	Capener Osteotome, Curved	260	12
PN0456	Capener Osteotome, Curved	260	15
PN0457	Capener Osteotome, Curved	260	20
PN0458	Capener Osteotome, Curved	260	25
PN0459	Capener Osteotome, Curved	260	30
PN0988	Capener Osteotome, Curved	260	40



PN0466

CODE	DESCRIPTION	LENGTH	WIDTH
PN0466	Capener Gouge, Straight	260	6
PN0986	Capener Gouge, Straight	260	8
PN0467	Capener Gouge, Straight	260	10
PN0044	Capener Gouge, Straight	260	12
PN0468	Capener Gouge, Straight	260	15
PN0469	Capener Gouge, Straight	260	20
PN0470	Capener Gouge, Straight	260	25
PN0699	Capener Gouge, Straight	260	30



PN0472

CODE	DESCRIPTION	LENGTH	WIDTH
PN0316	Capener Lamina Gouge, Cranked	240	6
PN0317	Capener Lamina Gouge, Cranked	240	8
PN1097	Capener Lamina Gouge, Cranked	240	10
PN0471	Capener Lamina Gouge, Cranked	240	11
PN1098	Capener Lamina Gouge, Cranked	240	13
PN0472	Capener Lamina Gouge, Cranked	240	16
PN1099	Capener Lamina Gouge, Cranked	240	19



PN0574

CODE	DESCRIPTION	LENGTH	WIDTH
PN1091	Swedish Osteotome	205	4
PN0816	Swedish Osteotome	205	6
PN0883	Swedish Osteotome	205	8
PN0574	Swedish Osteotome	205	10
PN0556	Swedish Osteotome	205	12
PN0527	Swedish Osteotome	205	15
PN0575	Swedish Osteotome	205	18
PN0524	Swedish Osteotome	205	20
PN0525	Swedish Osteotome	205	25
PN0526	Swedish Osteotome	205	30



PN0597

CODE	DESCRIPTION	LENGTH	WIDTH
PN0907	Swedish Chisel	205	10
PN0597	Swedish Chisel	205	12
PN0521	Swedish Chisel	205	15
PN0522	Swedish Chisel	205	20
PN0523	Swedish Chisel	205	25
PN1182	Swedish Chisel	205	30

SWEDISH GOUGE



PN0082

CODE	DESCRIPTION	LENGTH	WIDTH
PN0873	Swedish Gouge	205	10
PN0082	Swedish Gouge	205	12
PN0956	Swedish Gouge	205	15
PN0874	Swedish Gouge	205	20
PN0957	Swedish Gouge	205	25



PN0613



PN0602

CODE	DESCRIPTION	LENGTH	WIDTH
PN0607	Smith Petersen Osteotome, Straight	205	6
PN0608	Smith Petersen Osteotome, Straight	205	10
PN0609	Smith Petersen Osteotome, Straight	205	13
PN0610	Smith Petersen Osteotome, Straight	205	15
PN0611	Smith Petersen Osteotome, Straight	205	20
PN0612	Smith Petersen Osteotome, Straight	205	25
PN0613	Smith Petersen Osteotome, Straight	205	30
PN0600	Smith Petersen Osteotome, Curved	205	6
PN0601	Smith Petersen Osteotome, Curved	205	10
PN0602	Smith Petersen Osteotome, Curved	205	13
PN0603	Smith Petersen Osteotome, Curved	205	15
PN0604	Smith Petersen Osteotome, Curved	205	20
PN0605	Smith Petersen Osteotome, Curved	205	25
PN0606	Smith Petersen Osteotome, Curved	205	30



PN1055

CODE	DESCRIPTION	LENGTH	WIDTH
PN1051	Smith Petersen Gouge, Straight	205	6
PN1052	Smith Petersen Gouge, Straight	205	10
PN1053	Smith Petersen Gouge, Straight	205	13
PN1054	Smith Petersen Gouge, Straight	205	15
PN1055	Smith Petersen Gouge, Straight	205	20
PN1056	Smith Petersen Gouge, Straight	205	25
PN1057	Smith Petersen Gouge, Straight	205	30



PN0709

CODE	DESCRIPTION	LENGTH	WIDTH
PN1014	Hibbs Osteotome, Straight	204	4
PN1015	Hibbs Osteotome, Straight	204	5
PN1016	Hibbs Osteotome, Straight	204	6
PN1017	Hibbs Osteotome, Straight	204	7
PN1018	Hibbs Osteotome, Straight	204	8
PN1019	Hibbs Osteotome, Straight	204	9
PN1020	Hibbs Osteotome, Straight	204	10
PN1021	Hibbs Osteotome, Straight	204	11
PN1022	Hibbs Osteotome, Straight	204	12
PN1023	Hibbs Osteotome, Straight	204	13
PN1024	Hibbs Osteotome, Straight	204	14
PN1025	Hibbs Osteotome, Straight	204	15
PN1026	Hibbs Osteotome, Straight	204	16
PN1027	Hibbs Osteotome, Straight	204	17
PN0709	Hibbs Osteotome, Straight	204	18
PN1028	Hibbs Osteotome, Straight	204	19
PN1029	Hibbs Osteotome, Straight	204	20
PN1030	Hibbs Osteotome, Straight	204	25
PN1031	Hibbs Osteotome, Straight	204	30



PN0688

CODE	DESCRIPTION	LENGTH	WIDTH
PN1032	Hibbs Osteotome, Curved	203	4
PN1033	Hibbs Osteotome, Curved	203	5
PN1034	Hibbs Osteotome, Curved	203	6
PN1035	Hibbs Osteotome, Curved	203	7
PN1036	Hibbs Osteotome, Curved	203	8
PN1037	Hibbs Osteotome, Curved	203	9
PN1038	Hibbs Osteotome, Curved	203	10
PN1039	Hibbs Osteotome, Curved	203	11
PN1040	Hibbs Osteotome, Curved	203	12
PN1041	Hibbs Osteotome, Curved	203	13
PN1042	Hibbs Osteotome, Curved	203	14
PN1043	Hibbs Osteotome, Curved	203	15
PN1044	Hibbs Osteotome, Curved	203	16
PN1045	Hibbs Osteotome, Curved	203	17
PN1046	Hibbs Osteotome, Curved	203	18
PN1047	Hibbs Osteotome, Curved	203	19
PN1048	Hibbs Osteotome, Curved	203	20
PN1049	Hibbs Osteotome, Curved	203	25
PN0688	Hibbs Osteotome, Curved	203	30

LAMBOTTE MINI OSTEOTOME



PN0895

CODE	DESCRIPTION	LENGTH	WIDTH
PN0891	Lambotte Mini Osteotome, Straight	125	2
PN0892	Lambotte Mini Osteotome, Straight	125	4
PN0893	Lambotte Mini Osteotome, Straight	125	6
PN0894	Lambotte Mini Osteotome, Straight	125	8
PN0895	Lambotte Mini Osteotome, Straight	125	10

MINI OSTEOTOME



PN0900

CODE	DESCRIPTION	LENGTH	WIDTH
PN0896	Mini Osteotome	155	4
PN0897	Mini Osteotome	155	6
PN0898	Mini Osteotome	155	8
PN0899	Mini Osteotome	155	10
PN0900	Mini Osteotome	155	12



PN0817

CODE	DESCRIPTION	LENGTH	WIDTH
PN0904	East Grinstead Chisel	180	3
PN0817	East Grinstead Chisel	180	5

EAST GRINSTEAD OSTEOTOME



PN0818

CODE	DESCRIPTION	LENGTH	WIDTH
PN0905	East Grinstead Osteotome	180	3
PN0818	East Grinstead Osteotome	180	5



PN0728



PN0724



PN0725

CODE	DESCRIPTION	LENGTH	WIDTH
PN0728	Wagner Cement Chisel, Straight, Negative	335	9
PN0724	Wagner Cement Chisel, Angled Forward 10°, Negative	335	9
PN0725	Wagner Cement Chisel, Angled Backward 10°, Negative	335	9



PN0662



PN0663

CODE	DESCRIPTION	LENGTH	WIDTH
PN0662	Osteotome, 30°, Angled	297	16
PN0663	Osteotome, 30°, Angled, Crescentic	297	16



PN0786

CODE	DESCRIPTION	LENGTH	WIDTH
PN0960	Cobbs Elevator	280	10
PN0786	Cobbs Elevator	280	13
PN0787	Cobbs Elevator	280	19

ROWE ZYGOMATIC ELEVATOR



PN0233

CODE	DESCRIPTION	LENGTH
PN0233	Rowe Zygomatic Elevator	270

INGE LAMINA SPREADER



PN0695

CODE	DESCRIPTION	LENGTH	WIDTH
PN0695	Inge Lamina Spreader	175	6



PN0221



PN0222

CODE	DESCRIPTION	LENGTH	WIDTH
PN0221	Anterior Retractor	410	22
PN0222	Anterior Retractor	340	43



PN0226



PN0216



PN0217

CODE	DESCRIPTION	LENGTH	WIDTH
PN0226	Hohmann Retractor	230	18
PN0216	Hohmann Retractor	265	24
PN0217	Hohmann Retractor	240	43

LANGENBECK RETRACTOR



PN0223

CODE	DESCRIPTION	LENGTH	WIDTH
PN0223	Langenbeck Retractor	253	40

VOLKMANN SKIN RETRACTOR



PN0224

CODE	DESCRIPTION	LENGTH	WIDTH
PN0224	Volkmann Skin Retractor	205	36



PN0263



PN0229



PN0230

CODE	DESCRIPTION	LENGTH	WIDTH
PN0263	Müller Femoral Retractor	355	28
PN0229	Telescopic Arm for Müller Femoral Retractor	Open 660, Closed 365	
PN0230	Weight for Müller Femoral Retractor		

SHARP BONE HOOK



PN0231

CODE	DESCRIPTION	LENGTH
PN0231	Sharp Bone Hook	233

HOOK FOR CEMENT EXTRACTION



PN0220

CODE	DESCRIPTION	LENGTH
PN0220	Hook for Cement Extraction	343



PN0718



PN0719

CODE	DESCRIPTION	LENGTH
PN0718	Hook for Cement Extraction, Small	345
PN0719	Hook for Cement Extraction, Medium	345

OSTEOTOME



PN0245

CODE	DESCRIPTION	LENGTH	WIDTH
PN0245	Osteotome	280	16

CHISEL



PN0247

CODE	DESCRIPTION	LENGTH	WIDTH
PN0247	Chisel, Curved	283	16



PN0256

CODE	DESCRIPTION	LENGTH	WIDTH
PN0256	Swan Neck Gouge	285	25

CHISEL



PN0252

CODE	DESCRIPTION	LENGTH	WIDTH
PN0252	Chisel, Short Curved	280	9



PN0246



PN0248

CODE	DESCRIPTION	LENGTH	WIDTH
PN0246	Gouge, Double Angled	280	9
PN0248	Gouge, Straight	280	9



PN0262



PN0225



PN0473

CODE	DESCRIPTION	LENGTH
PN0262	Cup Introducer	360
PN0225	Antero-Lateral Cup Positioning Guide	
PN0473	Postero-Lateral Cup Positioning Guide	



1/1



1/1

CODE	DESCRIPTION	LENGTH	CUP SIZE
PN0259	Müller Curette, Small	270	9 x 13
PN0260	Müller Curette, Medium	275	13 x 20



PN0261



1/1



PN0258



1/1

CODE	DESCRIPTION	LENGTH	CUP SIZE
PN0261	Müller Curette, Large	240	27 x 37
PN0258	Müller Curette, Long	350	9 x 13

JOINT CAPSULE SCISSORS



PN0474

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0474	Joint Capsule Scissors	250	50

JOINT CAPSULE CLAMP



PN0228

CODE	DESCRIPTION	LENGTH	JAW LENGTH
PN0228	Joint Capsule Clamp	210	25



PN0215

CODE	DESCRIPTION	LENGTH	WIDTH
PN0215	Box Chisel	290	31

FEMORAL HEAD EXTRACTOR



PN0277

CODE	DESCRIPTION	LENGTH
PN0277	Femoral Head Extractor	230

FEMORAL DISLOCATION LEVER



PN0242

CODE	DESCRIPTION	LENGTH	WIDTH
PN0242	Femoral Dislocation Lever	285	26

FEMORAL STEM DRIVER



PN0251

CODE	DESCRIPTION	LENGTH	WIDTH
PN0251	Femoral Stem Driver	285	17



PN0283



PN0475

CODE	DESCRIPTION	LENGTH	SIZE
PN0283	Femoral Head Impactor	204	
PN0475	Nylon Head for Femoral Head Impactor		30 x 14

QUICK COUPLING HANDLE



PN0722

CODE	DESCRIPTION	LENGTH
PN0722	Quick Coupling Handle	191

CEMENT EXTRACTOR



PN0736

CODE	DESCRIPTION	SIZE
PN0736	Cement Extractor	8



PN0744

CODE	DESCRIPTION	SIZE
PN0737	Cement Reamer	7
PN0738	Cement Reamer	8
PN0739	Cement Reamer	9
PN0740	Cement Reamer	10
PN0741	Cement Reamer	11
PN0742	Cement Reamer	12
PN0743	Cement Reamer	13
PN0744	Cement Reamer	14

CENTERING DRILL



PN0735

CODE	DESCRIPTION	SIZE
PN0734	Centering Drill	6
PN0735	Centering Drill	8



PN0250



PN0287

CODE	DESCRIPTION	LENGTH	WIDTH
PN0250	Cement Chisel, Angled 15°	345	9
PN0287	Cement Splitting Chisel	350	11.5



PN0255



PN0254

CODE	DESCRIPTION	LENGTH	WIDTH
PN0255	Cement Chisel	325	5.5
PN0254	Cement Chisel	330	7.5



PN0253



PN0249



PN0257

CODE	DESCRIPTION	LENGTH	WIDTH
PN0253	Cement Chisel, Positive	345	9
PN0249	Cement Chisel, Positive	345	11
PN0257	Cement Chisel, Negative	265	9



PN0946



PN0944

CODE	DESCRIPTION	LENGTH	WIDTH
PN0946	Cement Chisel, Angled 15°, with Strike Plate	355	9
PN0944	Cement Splitting Chisel, with Strike Plate	360	11.5



PN0945



PN0732

CODE	DESCRIPTION	LENGTH	WIDTH
PN0945	Cement Chisel, with Strike Plate	335	5.5
PN0732	Cement Chisel, with Strike Plate	340	7.5



PN0731



PN0947

CODE	DESCRIPTION	LENGTH	WIDTH
PN0731	Cement Chisel, Positive, with Strike Plate	355	9
PN0947	Cement Chisel, Positive, with Strike Plate	355	11



PN0729

CODE	DESCRIPTION	LENGTH	WIDTH
PN0729	Cement Splitting Chisel	370	11.5

TROCHANTER FORCEPS



PN0227

CODE	DESCRIPTION	LENGTH
PN0227	Trochanter Forceps	215

RONGEUR FORCEPS



PN0218

CODE	DESCRIPTION	LENGTH
PN0218	Rongeur Forceps for Cement Extraction	278



PN0330

CODE	DESCRIPTION	LENGTH	WIDTH
PN0713	Osteotome, Thin Blade	150	6
PN0219	Osteotome, Thin Blade	150	8
PN0330	Osteotome, Thin Blade	150	10



PN0480

CODE	DESCRIPTION	LENGTH	WIDTH
PN0477	Chisel	305	5
PN0286	Chisel	305	10
PN0478	Chisel	305	15
PN0479	Chisel	305	20
PN0480	Chisel	305	25
PN0481	Chisel	305	30



PN0483

CODE	DESCRIPTION	LENGTH	WIDTH
PN0482	Osteotome	305	5
PN0483	Osteotome	305	10
PN0484	Osteotome	305	15
PN0485	Osteotome	305	20
PN0486	Osteotome	305	25
PN0244	Osteotome	305	30



PN0327

CODE	DESCRIPTION	LENGTH	WIDTH
PN0487	Gouge	305	5
PN0325	Gouge	305	10
PN0326	Gouge	305	15
PN0327	Gouge	305	20
PN0328	Gouge	305	25
PN0329	Gouge	305	30



PN0280



PN0281

CODE	DESCRIPTION	LENGTH
PN0280	Patella Resection Clamp	230
PN0281	Patella Impaction Clamp	230

WILKINSON ORTHOPAEDIC LOCKING FORCEPS



PN0292

CODE	DESCRIPTION	LENGTH
PN0292	Wilkinson Orthopaedic Locking Forceps	215

WILKINSON PUNCH



PN0214

CODE	DESCRIPTION	LENGTH
PN0214	Wilkinson Punch	240



PN0434



PN0398

CODE	DESCRIPTION	LENGTH	SIZE
PN0434	Charnley Hand Operated Femoral Reamer, Small	293	9.3 tapering to 6.3
PN0398	Charnley Hand Operated Femoral Reamer, Large	330	12.3 tapering to 9.5

CHARNLEY INITIAL INCISION RETRACTOR



PN0149

CODE	DESCRIPTION
PN0149	Charnley Initial Incision Retractor

Individual parts can be purchased separately:

CODE	DESCRIPTION	SIZE
PN1134	Charnley Initial Incision Retractor Frame	
PN0813	Charnley Initial Incision Retractor Weight & Chain	
PN1050	Charnley Initial Incision Retractor Blades, Standard	Length 62 x Depth 50 & Length 180 x Depth 40



PN0811




PN0812

CODE	DESCRIPTION	SIZE
PN0811	Charnley Initial Incision Retractor Blades, Long Reach	Length 62 x Depth 85 & Length 180 x Depth 75
PN1172	Charnley Initial Incision Retractor Blades, Extra Long Reach	Length 62 x Depth 135 & Length 180 x Depth 125
PN0812	Charnley Initial Incision Retractor Arm	

CHARNLEY RETRACTOR



PN0003



PN0596

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0003	Charnley Horizontal Retractor	285	18	Semi-Sharp
PN0596	Hardinge Horizontal Retractor	285	18	Semi-Sharp



PN0002



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0002	Charnley Vertical Retractor	245	22	3 x 4 Sharp



PN0264



1/1

PN0265



1/1

PN0264



1/1

PN0241

CODE	DESCRIPTION	LENGTH	CUP SIZE
PN0265	Charnley Bone Curette, T Handle, Small	370	11 x 20
PN0264	Charnley Bone Curette, T Handle, Medium	370	15 x 28
PN0241	Charnley Bone Curette, T Handle, Large	370	18 x 34



PN0408

CODE	DESCRIPTION	LENGTH
PN0408	Charnley Ring Curette	370

MORRIS RETRACTOR



PN0814

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0814	Morris Retractor	256	85	Sharp

OLLIER RETRACTOR



PN0815

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0815	Ollier Retractor	245	65	Blunt



PN0620



PN0614

CODE	DESCRIPTION	LENGTH	WIDTH
PN0620	Charnley Hohmann Bone Elevator, Straight	245	16
PN0614	Charnley Hohmann Bone Elevator, 45°	235	16

CHARNLEY BONE LEVER



PN0593

CODE	DESCRIPTION	LENGTH	WIDTH
PN0593	Charnley Bone Lever / Periosteal Elevator	210	10

CHARNLEY FEMORAL LEVER



PN0841

CODE	DESCRIPTION	LENGTH	WIDTH
PN0841	Charnley Femoral Lever	300	25



PN0121

CODE	DESCRIPTION	LENGTH	WIDTH
PN0121	Watson Jones Bone Lever	284	27

WATSON JONES GOUGE



PN0331

CODE	DESCRIPTION	LENGTH	WIDTH
PN0331	Watson Jones Gouge	305	15

FAGGE BONE LEVER



PN0120

CODE	DESCRIPTION	LENGTH	WIDTH
PN0120	Fagge Bone Lever	275	20

MURPHY SKID BONE LEVER



PN0165

CODE	DESCRIPTION	LENGTH	WIDTH
PN0165	Murphy Skid Bone Lever	300	30 and 20



PN0158



1/1



PN0692



PN0693

CODE	DESCRIPTION	COLOUR	LENGTH	JAW LENGTH
PN0158	Charnley Wire Holding Forceps	Plain	185	18
PN0692	Charnley Wire Holding Forceps	Blue	185	18
PN0693	Charnley Wire Holding Forceps	Black	185	18



PN0674



PN0675

CODE	DESCRIPTION	LENGTH
PN0674	Hardinge Cement Restrictor, Short	230
PN0675	Hardinge Cement Restrictor, Long	335



PN0723

CODE	DESCRIPTION	LENGTH
PN0723	Prosthesis Punch	255

WROBLEWSKI DISLOCATING HOOK



PN0745

CODE	DESCRIPTION	LENGTH
PN0745	Wroblewski Dislocating Hook	260



PN1063



PN1061

CODE	DESCRIPTION	LENGTH	SIZE
PN1113	Wroblewski Femoral Bone Block Kit - Punch	233	8
PN1062	Wroblewski Femoral Bone Block Kit - Punch	233	10
PN1063	Wroblewski Femoral Bone Block Kit - Punch	233	12
PN1064	Wroblewski Femoral Bone Block Kit - Punch	233	14
PN1061	Wroblewski Femoral Bone Block Kit - Trephine	152	



PN0543



PN0542

CODE	DESCRIPTION	LENGTH
PN0543	Hey Groves Bone Holding Forceps	210
PN0542	Hey Groves Bone Holding Forceps	310

Alphabetical Index

Description	Page	Description	Page
A		Charnley Femoral Lever	218
Anterior Retractor	175	Charnley Hand Operated Femoral Reamer	209
Antero-Lateral Cup Positioning Guide	185	Charnley Hohmann Bone Elevator	217
		Charnley Horizontal Retractor	212
B		Charnley Initial Incision Retractor	210-211
Beck Aortic Clamp	70	Charnley Ring Curette	215
Bone Block Kit - Wroblewski Femoral	224	Charnley Vertical Retractor	213
Bone Cutting Forceps - Horsley Liston	156	Charnley Wire Holding Forceps	221
Bone Cutting Forceps - Liston	150-151	Chisel - Box	189
Bone Cutting Forceps - McIndoe	149	Chisel - Capener	158
Bone Cutting Forceps - Paton	148	Chisel - Cement	194-201
Bone Cutting Forceps - Stamm	148	Chisel - Curved 16mm	182
Bone Cutting Forceps - Ward	153	Chisel - East Grinstead	169
Bone Nibbling Forceps - Cawthorne	154	Chisel - Extra Long	204
Box Chisel	189	Chisel - Short Curved, 9mm	183
Brock Auricular Clamp	86	Chisel - Swedish	163
Browse Retractor	22	Chisel - Wagner	171
Bulldog Clamp - Cooley, Cross Action	32-33	Clamp - Beck Aortic	70
Bulldog Clamp - Cooley, Ring Handle	40	Clamp - Brock Auricular	86
Bulldog Clamp - Debakey, Cross Action	34-35	Clamp - Castenada	52
Bulldog Clamp - Debakey, Ring Handle	41	Clamp - Cooley Bulldog, Cross Action	32-33
Bulldog Clamp - Wolfe Modified Cooley, Ring Handle	42	Clamp - Cooley Bulldog, Ring Handle	40
		Clamp - Cooley Derra	38
C		Clamp - Cooley Multi Purpose	44
Cairns Scalp Retractor	24	Clamp - Cooley Paediatric Occlusion	36
Capener Chisel	158	Clamp - Cooley Vessel	42
Capener Gouge, Lamina	161	Clamp - Crafoord Coarctation	82
Capener Gouge, Straight	160	Clamp - Crafoord Dissecting	83
Capener Osteotome	159	Clamp - Dardik Multi Purpose	50-51
Castenada Clamp	52	Clamp - Debakey Aneurysm	60-61
Cawthorne Bone Nibbling Forceps	154	Clamp - Debakey Aortic	66-67
Cawthorne Retractor	12	Clamp - Debakey Aortic Aneurysm	62-64
Cement Chisel	194-201	Clamp - Debakey Aortic Occlusion	65
Cement Extraction Hook	180-181	Clamp - Debakey Bainbridge	75-76
Cement Extraction Rongeur Forceps	202	Clamp - Debakey Bulldog, Cross Action	34-35
Cement Extractor	192	Clamp - Debakey Bulldog, Ring Handle	41
Cement Reamer	193	Clamp - Debakey Coarctation	77-80
Cement Restrictor - Hardinge	222	Clamp - Debakey Derra	39
Centering Drill	193	Clamp - Debakey Hysterectomy	124-127
Charnley Bone Curette	214	Clamp - Debakey Iliac	57
Charnley Bone Lever / Periosteal Elevator	218	Clamp - Debakey Miniature	43

Alphabetical Index

Description	Page	Description	Page
Clamp - Debakey Multi Purpose	45-49	Clamp - Wolfe Modified Cooley Bulldog, Ring Handle	42
Clamp - Debakey Paediatric Occlusion	37	Cobbs Elevator	173
Clamp - Debakey Peripheral	54-55	Collingwood Stewart Hernia Forceps	109
Clamp - Debakey Renal	57	Cone Retractor	18-20
Clamp - Debakey Tangential	58-59	Cooley Bulldog Clamp, Cross Action	32-33
Clamp - Doyen	112-115	Cooley Bulldog Clamp, Ring Handle	40
Clamp - Einsell Foss Colon	135	Cooley Derra Clamp	38
Clamp - Finch Excision	134	Cooley Multi Purpose Clamp	44
Clamp - Fistula	111	Cooley Paediatric Occlusion Clamp	36
Clamp - Glover	84-85	Cooley Vessel Clamp	42
Clamp - Goldberg Anterior Resection	137	Cottle Nasal Speculum	142-143
Clamp - Gwilliams Hysterectomy	130	Crafoord Coarctation Clamp	82
Clamp - Hayes Colon Resection	136	Crafoord Dissecting Clamp	83
Clamp - Herrick Kidney Pedicle	111	Cup Introducer	185
Clamp - Houghton Excision	135	Cup Positioning Guide - Antero-Lateral	185
Clamp - Howkin Hysterectomy	133	Cup Positioning Guide - Postero-Lateral	185
Clamp - Javid Carotid Artery	74	Curette - Charnley Bone	214
Clamp - Joint Capsule	188	Curette - Charnley Ring	215
Clamp - Lambert Kay	69	Curette - Müller	186-187
Clamp - Lane Twin	120-121		
Clamp - Lees Bronchus	87	D	
Clamp - Leland Jones Peripheral Vascular	56	Dardik Multi Purpose Clamp	50-51
Clamp - Mattox Aorta	65	Debakey Aneurysm Clamp	60-61
Clamp - McDougall Prostatectomy	138	Debakey Aortic Aneurysm Clamp	62-64
Clamp - Morris Aorta	68	Debakey Aortic Clamp	66-67
Clamp - Moynihan	102	Debakey Aortic Occlusion Clamp	65
Clamp - Neonatal Miniature	53	Debakey Bainbridge Clamp	75-76
Clamp - Parker Kerr	116-117	Debakey Bulldog Clamp, Cross Action	34-35
Clamp - Patella	207	Debakey Bulldog Clamp, Ring Handle	41
Clamp - Payr Intestinal	118	Debakey Coarctation Clamp	77-80
Clamp - Potts Coarctation	81	Debakey Derra Clamp	39
Clamp - Price Thomas Bronchus	88	Debakey Hysterectomy Clamp	124-127
Clamp - Read Stallworthy Hysterectomy	132	Debakey Iliac Clamp	57
Clamp - Rogers Hysterectomy	131	Debakey Ligature Forceps	91
Clamp - Ross Vena Cava	73	Debakey Miniature Clamp	43
Clamp - Satinsky	71-72	Debakey Modified Semb Ligature Passing Forceps	91
Clamp - Swinney	94	Debakey Multi Purpose Clamp	45-49
Clamp - Tanner (Lane) Twin	119	Debakey Paediatric Occlusion Clamp	37
Clamp - Tubing	95-97	Debakey Peripheral Clamp	54-55
Clamp - Ulrich	92-93	Debakey Renal Clamp	57
Clamp - Waterston Vessel	108	Debakey Tangential Clamp	58-59

Alphabetical Index

Description	Page	Description	Page
Doyen Clamp	112-115	Forceps - Nelson Roberts	101
Drill - Centering	193	Forceps - Purse String	110
E		Forceps - Roberts Artery	102
East Grinstead Chisel	169	Forceps - Sawtell	101
East Grinstead Osteotome	169	Forceps - Semb Ligature Passing	89-90
Einsell Foss Colon Clamp	135	Forceps - Teale Vulsellum	133
Eisenhammer Retractor	26	Forceps - Trochanter	202
Elevator - Charnley Bone Lever / Periosteal	218	Forceps - Waterston Dissecting	103-107
Elevator - Charnley Hohmann Bone	217	Forceps - Wilkinson Orthopaedic Locking	208
Elevator - Cobbs	173	G	
Elevator - Hohmann	176-177	Geary Grant Forceps	122-123
Elevator - Rowe Zygomatic	174	Gelpi Retractor	14-15
Extractor - Cement	192	Glasgow Rongeur	152
Extractor - Femoral Head	189	Glover Clamp	84-85
F		Goldberg Anterior Resection Clamp	137
Fagge Bone Lever	220	Gouge - Capener, Lamina	161
Farrands Retractor	27	Gouge - Capener, Straight	160
Femoral Dislocation Lever	190	Gouge - Double Angled, 9mm	184
Femoral Head Extractor	189	Gouge - Extra Long	206
Femoral Head Impactor	191	Gouge - Smith Petersen	165
Femoral Lever - Charnley	218	Gouge - Straight, 9mm	184
Femoral Stem Driver	190	Gouge - Swan Neck	183
Finch Excision Clamp	134	Gouge - Swedish	163
Fistula Clamp	111	Gouge - Watson Jones	219
Forceps - Cement Extraction Rongeur	202	Gwilliams Hysterectomy Clamp	130
Forceps - Charnley Wire Holding	221	H	
Forceps - Collingwood Stewart Hernia	109	Handle - Quick Coupling	192
Forceps - Debakey Ligature	91	Hardinge Cement Restrictor	222
Forceps - Debakey Modified Semb Ligature Passing	91	Hardinge Horizontal Retractor	212
Forceps - Geary Grant	122-123	Hardman Retractor	16-17
Forceps - Hey Groves Bone Holding	225	Harefield Sternal Blades	29
Forceps - Irwin Moore Nasal	146	Harvey Jackson Laminectomy Retractor	21
Forceps - King Vein Graft Press	110	Hayes Colon Resection Clamp	136
Forceps - Lahey Artery	98-99	Henderson Retractor	8
Forceps - Lincoln Dissecting	109	Herrick Kidney Pedicle Clamp	111
Forceps - Lloyd Davis Artery	100	Hey Groves Bone Holding Forceps	225
Forceps - Lloyd Davis Rectal Occlusion	134	Hibbs Osteotome	166-167
Forceps - Luc Nasal	144-145	Hohmann Retractor	176-177
Forceps - Maingot Hysterectomy	128-129	Hook - Cement Extraction	180-181

Alphabetical Index

Description	Page	Description	Page
Hook - Phlebectomy	30	Lees Bronchus Clamp	87
Hook - Sharp Bone	180	Leland Jones Peripheral Vascular Clamp	56
Hook - Wroblewski Dislocating	223	Lever - Charnley Bone / Periosteal Elevator	218
Horsley Liston Bone Cutting Forceps	156	Lever - Charnley Femoral	218
Houghton Excision Clamp	135	Lever - Fagge Bone	220
Howkin Hysterectomy Clamp	133	Lever - Femoral Dislocation	190
Hysterectomy - Debakey	124-127	Lever - Murphy Skid Bone	220
Hysterectomy - Gwilliams	130	Lever - Watson Jones Bone	219
Hysterectomy - Howkin	133	Lincoln Dissecting Forceps	109
Hysterectomy - Maingot	128-129	Liston Bone Cutting Forceps	150-151
Hysterectomy - Read Stallworthy	132	Lloyd Davis Artery Forceps	100
Hysterectomy - Rogers	131	Lloyd Davis Rectal Occlusion Forceps	134
		Luc Nasal Forceps	144-145
		Luer Jansen Rongeur	154
I		M	
IMA Retractor	30	Maingot Hysterectomy Forceps	128-129
Impactor - Femoral Head	191	Mattox Aorta Clamp	65
Inge Lamina Spreader	174	Mawson Retractor	12
Irwin Moore Nasal Forceps	146	Mayo Adson Laminectomy Retractor	21
Irwin Retractor	9	McDougall Prostatectomy Clamp	138
		McIndoe Bone Cutting Forceps	149
J		Mini Osteotome	168
Jackson Burrows Retractor	13	Mollison Retractor	5
Jansen Middleton Septum Forceps	155	Morris Aorta Clamp	68
Javid Carotid Artery Clamp	74	Morris Retractor	216
Jefferson Retractor	9	Moynihan Clamp	102
Joint Capsule Clamp	188	Müller Curette	186-187
Joint Capsule Scissors	188	Müller Femoral Retractor	179
Jolls Thyroid Retractor	7	Multi Purpose - Cooley	44
Juniper Retractor	13	Multi Purpose - Dardik	50-51
		Multi Purpose - Debakey	45-49
K		Murphy Skid Bone Lever	220
Killian Nasal Speculum	140-141		
King Vein Graft Press Forceps	110		
		N	
L		Nasal Forceps - Irwin Moore	146
Lahey Artery Forceps	98-99	Nasal Forceps - Luc	144-145
Lambert Kay Clamp	69	Nasal Speculum - Cottle	142-143
Lambotte Mini Osteotome	168	Nasal Speculum - Killian	140-141
Lamina Spreader - Inge	174	Nelson Roberts Forceps	101
Lane Twin Clamp	120-121	Neonatal Miniature Clamp	53
Langenbeck Retractor	178	Norfolk & Norwich Retractor	6-7

Alphabetical Index



Description	Page	Description	Page
O		Retractor - Cairns Scalp	24
Olivecrona Rongeur	155	Retractor - Cawthorne	12
Ollier Retractor	216	Retractor - Charnley Horizontal	212
Osteotome - 16mm	182	Retractor - Charnley Initial Incision	210-211
Osteotome - 30°, Angled	172	Retractor - Charnley Vertical	213
Osteotome - 30°, Angled, Crescentic	172	Retractor - Cone	18-20
Osteotome - Capener	159	Retractor - Eisenhammer	26
Osteotome - East Grinstead	169	Retractor - Farrands	27
Osteotome - Extra Long	205	Retractor - Gelpi	14-15
Osteotome - Hibbs	166-167	Retractor - Hardinge Horizontal	212
Osteotome - Lambotte Mini	168	Retractor - Hardman	16-17
Osteotome - Mini	168	Retractor - Harvey Jackson Laminectomy	21
Osteotome - Smith Petersen	164	Retractor - Henderson	8
Osteotome - Swedish	162	Retractor - Hohmann	176-177
Osteotome - Thin Blade	203	Retractor - IMA	30
		Retractor - Irwin	9
P		Retractor - Jackson Burrows	13
Park Anal Retractor	25	Retractor - Jefferson	9
Parker Kerr Clamp	116-117	Retractor - Jolls Thyroid	7
Patella Clamp	207	Retractor - Juniper	13
Paton Bone Cutting Forceps	148	Retractor - Langenbeck	178
Payr Intestinal Clamp	118	Retractor - Mawson	12
Pennybacker Rongeur	152	Retractor - Mayo Adson Laminectomy	21
Phlebectomy Hook	30	Retractor - Mollison	5
Plester Jansen Retractor	10-11	Retractor - Morris	216
Postero-Lateral Cup Positioning Guide	185	Retractor - Müller Femoral	179
Potts Coarctation Clamp	81	Retractor - Norfolk & Norwich	6-7
Price Thomas Bronchus Clamp	88	Retractor - Ollier	216
Prosthesis Punch	223	Retractor - Park Anal	25
Purse String Forceps	110	Retractor - Plester Jansen	10-11
		Retractor - R.N.T.N.E.H.	24
Q		Retractor - Sellors	28
Quick Coupling Handle	192	Retractor - Talbot Orbital	23
		Retractor - Travers	4
R		Retractor - Volkmann Skin	178
R.N.T.N.E.H. Retractor	24	Retractor - Wallace Modified Plester Jansen	11
Read Stallworthy Hysterectomy Clamp	132	Retractor - Weislander	2
Reamer - Cement	193	Retractor - Weitlaner	3
Reamer - Charnley Hand Operated Femoral	209	Retractor - West	2
Retractor - Anterior	175	Roberts Artery Forceps	102
Retractor - Browse	22	Rogers Hysterectomy Clamp	131

Alphabetical Index

Description	Page	Description	Page
Rongeur - Forceps for Cement Extraction	202	W	
Rongeur - Glasgow	152	Wagner Cement Chisel	170-171
Rongeur - Luer Jansen	154	Wallace Modified Plester Jansen Retractor	11
Rongeur - Olivecrona	155	Ward Bone Cutting Forceps	153
Rongeur - Pennybacker	152	Waterston Dissecting Forceps	103-107
Rongeur - Van Havre	156	Waterston Vessel Clamp	108
Rongeur - Wilms	153	Watson Jones Bone Lever	219
Ross Vena Cava Clamp	73	Watson Jones Gouge	219
Rowe Zygomatic Elevator	174	Weislander Retractor	2
		Weitlaner Retractor	3
S		West Retractor	2
Satinsky Clamp	71-72	Wilkinson Orthopaedic Locking Forceps	208
Sawtell Forceps	101	Wilkinson Punch	208
Scissors - Joint Capsule	188	Wilms Rongeur	153
Sellers Retractor	28	Wolfe Modified Cooley Bulldog Clamp, Ring Handle	42
Semb Ligature Passing Forceps	89-90	Wroblewski Dislocating Hook	223
Septum Forceps - Jansen Middleton	155	Wroblewski Femoral Bone Block Kit	224
Sharp Bone Hook	180		
Smith Petersen Gouge	165		
Smith Petersen Osteotome	164		
Stamm Bone Cutting Forceps	148		
Swan Neck Gouge	183		
Swedish Chisel	163		
Swedish Gouge	163		
Swedish Osteotome	162		
Swinney Clamp	94		
T			
Talbot Orbital Retractor	23		
Tanner (Lane) Twin Clamp	119		
Teale Vulsellum Forceps	133		
Travers Retractor	4		
Trochanter Forceps	202		
Tubing Clamp	95-97		
U			
Ulrich Clamp	92-93		
V			
Van Havre Rongeur	156		
Volkman Skin Retractor	178		

Numerical Index

Code	Page	Code	Page	Code	Page	Code	Page
PN0001	2	PN0045	46	PN0108	136	PN0167	144
PN0002	213	PN0046	91	PN0109	136	PN0168	144
PN0003	212	PN0052	80	PN0110	14	PN0169	144
PN0004	13	PN0055	137	PN0115	24	PN0170	145
PN0005	14	PN0057	35	PN0116	51	PN0171	145
PN0006	21	PN0059	53	PN0117	51	PN0172	145
PN0007	6	PN0061	134	PN0118	32	PN0174	99
PN0008	9	PN0062	53	PN0119	51	PN0175	99
PN0009	5	PN0063	53	PN0120	220	PN0176	40
PN0010	21	PN0064	56	PN0121	219	PN0177	52
PN0011	24	PN0067	15	PN0122	50	PN0178	40
PN0012	5	PN0068	70	PN0123	119	PN0179	40
PN0013	8	PN0069	72	PN0124	50	PN0180	65
PN0014	2	PN0070	72	PN0125	16	PN0181	36
PN0015	4	PN0071	72	PN0128	56	PN0182	36
PN0016	4	PN0077	78	PN0129	114	PN0183	36
PN0017	12	PN0078	78	PN0130	115	PN0184	36
PN0018	9	PN0079	84	PN0131	27	PN0185	36
PN0019	12	PN0082	163	PN0134	37	PN0188	41
PN0020	23	PN0083	48	PN0137	158	PN0189	41
PN0021	6	PN0084	132	PN0138	131	PN0191	102
PN0022	17	PN0085	66-67	PN0139	130	PN0211	22
PN0023	16	PN0086	89	PN0140	71	PN0212	7
PN0024	65	PN0087	110	PN0141	132	PN0213	118
PN0025	54-55	PN0088	86	PN0142	124	PN0214	208
PN0026	43	PN0089	159	PN0143	42	PN0215	189
PN0027	43	PN0092	142-143	PN0144	42	PN0216	176-177
PN0028	137	PN0094	151	PN0146	126	PN0217	177
PN0029	39	PN0095	109	PN0149	210	PN0218	202
PN0030	39	PN0096	85	PN0152	102	PN0219	203
PN0031	39	PN0097	81	PN0153	128	PN0220	180
PN0032	131	PN0098	103	PN0154	46	PN0221	175
PN0033	54-55	PN0099	61	PN0157	82	PN0222	175
PN0034	44	PN0100	103	PN0158	221	PN0223	178
PN0035	44	PN0101	60-61	PN0159	96	PN0224	178
PN0036	44	PN0102	62	PN0160	96	PN0225	185
PN0038	15	PN0103	64	PN0161	96	PN0226	176-177
PN0039	81	PN0104	109	PN0163	79	PN0227	202
PN0040	142-143	PN0105	67	PN0164	69	PN0228	188
PN0041	28	PN0106	66-67	PN0165	220	PN0229	179
PN0044	160	PN0107	63	PN0166	73	PN0230	179

Numerical Index

Code	Page	Code	Page	Code	Page	Code	Page
PN0231	180	PN0280	207	PN0331	219	PN0388	115
PN0233	174	PN0281	207	PN0332	150	PN0389	135
PN0236	150	PN0283	191	PN0333	152	PN0390	134
PN0237	155	PN0286	204	PN0335	111	PN0393	84
PN0238	156	PN0287	194	PN0336	38	PN0394	85
PN0241	214	PN0290	142-143	PN0337	45	PN0395	130
PN0242	190	PN0291	74	PN0338	149	PN0396	135
PN0243	88	PN0292	208	PN0339	86	PN0398	209
PN0244	205	PN0293	89	PN0340	86	PN0399	121
PN0245	182	PN0294	57	PN0341	40	PN0400	120
PN0246	184	PN0295	83	PN0342	40	PN0401	87
PN0247	182	PN0297	138	PN0347	82	PN0402	87
PN0248	184	PN0298	138	PN0348	83	PN0403	100
PN0249	196-197	PN0299	117	PN0349	50	PN0405	128
PN0250	194	PN0300	116	PN0350	50	PN0406	129
PN0251	190	PN0301	68	PN0351	50	PN0407	68
PN0252	183	PN0302	108	PN0352	51	PN0408	215
PN0253	196-197	PN0305	30	PN0353	51	PN0410	11
PN0254	195	PN0306	112	PN0355	49	PN0414	118
PN0255	195	PN0307	123	PN0356	49	PN0417	91
PN0256	183	PN0308	35	PN0357	133	PN0418	90
PN0257	197	PN0309	34	PN0358	62	PN0419	94
PN0258	187	PN0310	33	PN0359	76	PN0420	94
PN0259	186	PN0311	159	PN0360	77	PN0421	97
PN0260	186	PN0312	159	PN0361	77	PN0422	97
PN0261	187	PN0313	159	PN0364	125	PN0423	97
PN0262	185	PN0314	159	PN0365	124	PN0424	93
PN0263	179	PN0315	159	PN0366	125	PN0425	92-93
PN0264	214	PN0316	161	PN0367	127	PN0426	92-93
PN0265	214	PN0317	161	PN0368	126	PN0427	140-141
PN0267	95	PN0318	38	PN0369	127	PN0428	140-141
PN0268	146	PN0319	129	PN0370	43	PN0429	140-141
PN0270	117	PN0320	142-143	PN0371	45	PN0430	140-141
PN0272	38	PN0322	100	PN0372	47	PN0431	140-141
PN0273	88	PN0324	25	PN0373	47	PN0432	148
PN0274	88	PN0325	206	PN0375	41	PN0434	209
PN0275	152	PN0326	206	PN0376	41	PN0435	156
PN0276	154	PN0327	206	PN0382	55	PN0436	98
PN0277	189	PN0328	206	PN0383	111	PN0437	151
PN0278	149	PN0329	206	PN0384	60-61	PN0440	22
PN0279	148	PN0330	203	PN0387	114	PN0441	23

Numerical Index

Code	Page	Code	Page	Code	Page	Code	Page
PN0443	10	PN0485	205	PN0585	107	PN0695	174
PN0444	10	PN0486	205	PN0587	57	PN0696	104
PN0445	11	PN0487	206	PN0589	153	PN0699	160
PN0446	20	PN0492	51	PN0590	34	PN0702	159
PN0447	18	PN0494	53	PN0591	34	PN0705	81
PN0448	18	PN0495	20	PN0592	32	PN0706	105
PN0449	19	PN0507	79	PN0593	218	PN0709	166
PN0450	19	PN0521	163	PN0594	33	PN0710	70
PN0451	159	PN0522	163	PN0596	212	PN0711	104
PN0452	159	PN0523	163	PN0597	163	PN0713	203
PN0453	159	PN0524	162	PN0598	52	PN0718	181
PN0454	159	PN0525	162	PN0599	35	PN0719	181
PN0455	159	PN0526	162	PN0600	164	PN0722	192
PN0456	159	PN0527	162	PN0601	164	PN0723	223
PN0457	159	PN0532	74	PN0602	164	PN0724	170-171
PN0458	159	PN0540	48	PN0603	164	PN0725	171
PN0459	159	PN0541	110	PN0604	164	PN0728	170-171
PN0460	158	PN0542	225	PN0605	164	PN0729	201
PN0461	158	PN0543	225	PN0606	164	PN0731	200
PN0462	158	PN0544	154	PN0607	164	PN0732	199
PN0463	158	PN0545	34	PN0608	164	PN0734	193
PN0464	158	PN0546	35	PN0609	164	PN0735	193
PN0465	158	PN0547	15	PN0610	164	PN0736	192
PN0466	160	PN0549	17	PN0611	164	PN0737	193
PN0467	160	PN0550	37	PN0612	164	PN0738	193
PN0468	160	PN0551	37	PN0613	164	PN0739	193
PN0469	160	PN0552	116	PN0614	217	PN0740	193
PN0470	160	PN0554	35	PN0619	153	PN0741	193
PN0471	161	PN0555	56	PN0620	217	PN0742	193
PN0472	161	PN0556	162	PN0627	113	PN0743	193
PN0473	185	PN0557	13	PN0628	112	PN0744	193
PN0474	188	PN0558	26	PN0647	73	PN0745	223
PN0475	191	PN0564	133	PN0662	172	PN0746	69
PN0477	204	PN0569	41	PN0663	172	PN0786	173
PN0478	204	PN0571	81	PN0674	222	PN0787	173
PN0479	204	PN0573	26	PN0675	222	PN0801	30
PN0480	204	PN0574	162	PN0679	113	PN0802	30
PN0481	204	PN0575	162	PN0688	167	PN0803	30
PN0482	205	PN0578	107	PN0692	221	PN0811	211
PN0483	205	PN0583	108	PN0693	221	PN0812	211
PN0484	205	PN0584	73	PN0694	105	PN0813	210

Numerical Index

Code	Page	Code	Page	Code	Page	Code	Page
PN0814	216	PN0980	158	PN1045	167	PN1092	37
PN0815	216	PN0981	158	PN1046	167	PN1093	34
PN0816	162	PN0982	158	PN1047	167	PN1094	3
PN0817	169	PN0983	158	PN1048	167	PN1095	122
PN0818	169	PN0984	158	PN1049	167	PN1096	56
PN0841	218	PN0985	158	PN1050	210	PN1097	161
PN0873	163	PN0986	160	PN1051	165	PN1098	161
PN0874	163	PN0987	159	PN1052	165	PN1099	161
PN0879	155	PN0988	159	PN1053	165	PN1113	224
PN0883	162	PN0989	8	PN1054	165	PN1134	210
PN0891	168	PN1014	166	PN1055	165	PN1136	142-143
PN0892	168	PN1015	166	PN1056	165	PN1137	142-143
PN0893	168	PN1016	166	PN1057	165	PN1138	142-143
PN0894	168	PN1017	166	PN1058	29	PN1139	142-143
PN0895	168	PN1018	166	PN1059	29	PN1172	211
PN0896	168	PN1019	166	PN1061	224	PN1182	163
PN0897	168	PN1020	166	PN1062	224		
PN0898	168	PN1021	166	PN1063	224		
PN0899	168	PN1022	166	PN1064	224		
PN0900	168	PN1023	166	PN1069	52		
PN0904	169	PN1024	166	PN1070	52		
PN0905	169	PN1025	166	PN1071	52		
PN0907	163	PN1026	166	PN1072	52		
PN0918	7	PN1027	166	PN1073	52		
PN0944	198	PN1028	166	PN1074	52		
PN0945	199	PN1029	166	PN1075	52		
PN0946	198	PN1030	166	PN1076	52		
PN0947	200	PN1031	166	PN1077	53		
PN0948	106	PN1032	167	PN1078	58-59		
PN0949	106	PN1033	167	PN1079	58-59		
PN0951	3	PN1034	167	PN1080	58-59		
PN0952	3	PN1035	167	PN1081	58-59		
PN0956	163	PN1036	167	PN1082	58-59		
PN0957	163	PN1037	167	PN1084	32		
PN0960	173	PN1038	167	PN1085	32		
PN0965	101	PN1039	167	PN1086	33		
PN0972	50	PN1040	167	PN1087	33		
PN0973	101	PN1041	167	PN1088	75		
PN0977	158	PN1042	167	PN1089	76		
PN0978	158	PN1043	167	PN1090	75		
PN0979	158	PN1044	167	PN1091	162		





Platts & Nisbett Ltd Woodfold Works Woodfold Sheffield S3 9PE
T: 0114 275 0387 **F:** 0114 279 8434 **E:** sales@plattsnisbett.com **W:** www.plattsnisbett.com
Registered in England: 1297203